


SN-001

Chemical Properties

CAS No. : 727699-84-5
 Formula: C₂₆H₂₁FN₂O₄S
 Molecular Weight: 476.52
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	SN-001 was reported as a STING inhibitor, which inhibited the activation of the STING signal pathway and to prevent or treat a STING- mediated disease.
Targets(IC50)	STING

Solubility Information

Solubility	DMSO: 60 mg/mL (125.91 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+90% Corn Oil: 2.5 mg/mL (5.25 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0985 mL	10.4927 mL	20.9855 mL
5 mM	0.4197 mL	2.0985 mL	4.1971 mL
10 mM	0.2099 mL	1.0493 mL	2.0985 mL
50 mM	0.042 mL	0.2099 mL	0.4197 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481