


Fadaltran

Chemical Properties

CAS No. :	1799809-36-1
Formula:	C ₂₄ H ₂₉ N ₅ O ₂
Molecular Weight:	419.52
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fadaltran is an α 2-adrenoreceptor antagonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3837 mL	11.9184 mL	23.8368 mL
5 mM	0.4767 mL	2.3837 mL	4.7674 mL
10 mM	0.2384 mL	1.1918 mL	2.3837 mL
50 mM	0.0477 mL	0.2384 mL	0.4767 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481