


Cbz-N-PEG10-acid

Chemical Properties

CAS No. :

Formula: C31H53NO14

Molecular Weight: 663.75



Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cbz-N-PEG10-acid is a PEG-based linker for PROTACs that connects two essential ligands, facilitating the formation of PROTAC molecules. This linker enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands connected by a linker: one ligand binds to an E3 ubiquitin ligase and the other to the target protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5066 mL	7.533 mL	15.0659 mL
5 mM	0.3013 mL	1.5066 mL	3.0132 mL
10 mM	0.1507 mL	0.7533 mL	1.5066 mL
50 mM	0.0301 mL	0.1507 mL	0.3013 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481