


Kartogenin

Chemical Properties

CAS No. :	4727-31-5
Formula:	C ₂₀ H ₁₅ NO ₃
Molecular Weight:	317.34
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Kartogenin (KGN) is an activator of the smad4/smad5 pathway, and promotes the selective differentiation of multipotent mesenchymal stem cells into chondrocytes.
Targets(IC50)	TGF-beta/Smad
In vitro	Kartogenin promotes type I collagen synthesis in dermal fibroblasts by activating the Smad4/Smad5 pathway, without exhibiting significant cytotoxicity on the activity, morphology, and survival of fibroblasts. Additionally, by enhancing the nuclear localization of CBFβ, Kartogenin induces differentiation of hMSCs into chondrocytes and protects articular chondrocytes.
In vivo	Kartogenin at a concentration of 100 nM stimulates collagen synthesis in the dermis of mice. Additionally, it promotes the growth of skeletal limbs and the development of synovial joints in mice. Concentrations of 10 μM and 100 μM Kartogenin respectively enhance cartilage repair in mouse models of osteoarthritis (OA) induced by collagenase VII and surgically.
Cell Research	For the cell viability assay, fibroblasts are seeded on six-well plates, followed by Kartogenin or DMSO treatment. According to the manufacturer's instructions, cell viability and viability are assessed by AlamarBlue assay.(Only for Reference)

Solubility Information

Solubility	H ₂ O: <1 mg/mL, DMSO: 41.67 mg/mL (131.31 mM),Sonication is recommended. Ethanol: 18 mg/mL (56.72 mM),Heating is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (6.3 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1512 mL	15.756 mL	31.5119 mL
5 mM	0.6302 mL	3.1512 mL	6.3024 mL
10 mM	0.3151 mL	1.5756 mL	3.1512 mL
50 mM	0.063 mL	0.3151 mL	0.6302 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Johnson K, et al. *Science*. 2012, 336(6082), 717-721.

Wang J, et al. *Biochem Biophys Res Commun*. 2014, 450(1), 568-574.

Decker RS, et al. *Dev Biol*. 2014, 395(2), 255-267.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481