


PI3K $\delta$ -IN-11

## Chemical Properties

CAS No. :	2413257-51-7
Formula:	C <sub>27</sub> H <sub>21</sub> N <sub>5</sub> O
Molecular Weight:	431.49
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	PI3K $\delta$ -IN-11 is a potent and selective PI3K $\delta$ inhibitor (IC <sub>50</sub> : 27.5 nM) that dose-dependently blocks PI3K/Akt pathway activity. PI3K $\delta$ -IN-11 can be used in B- or T-cell-associated malignancy studies.
Targets(IC <sub>50</sub> )	Apoptosis,Others,Akt,PI3K

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3176 mL	11.5878 mL	23.1755 mL
5 mM	0.4635 mL	2.3176 mL	4.6351 mL
10 mM	0.2318 mL	1.1588 mL	2.3176 mL
50 mM	0.0464 mL	0.2318 mL	0.4635 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481