


Edoxaban Isomer

Chemical Properties

CAS No. :	1255529-24-8
Formula:	C ₂₄ H ₃₀ ClN ₇ O ₄ S
Molecular Weight:	548.06
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Edoxaban Isomer is an impurity of Edoxaban, which is a novel inhibitor of factor Xa.
Targets(IC50)	Factor Xa,Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	-----------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8246 mL	9.1231 mL	18.2462 mL
5 mM	0.3649 mL	1.8246 mL	3.6492 mL
10 mM	0.1825 mL	0.9123 mL	1.8246 mL
50 mM	0.0365 mL	0.1825 mL	0.3649 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481