


HPK1-IN-26

## Chemical Properties

CAS No. : 2229042-24-2  
 Formula: C<sub>19</sub>H<sub>21</sub>N<sub>5</sub>O<sub>5</sub>  
 Molecular Weight: 367.47  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	HPK1-IN-26, a novel compound, is a potent inhibitor of both HPK1 and GLK. This small molecule exhibits promise in aiding the study of animal pathogen infection [1].
Targets(IC50)	Others,MAPK

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7213 mL	13.6066 mL	27.2131 mL
5 mM	0.5443 mL	2.7213 mL	5.4426 mL
10 mM	0.2721 mL	1.3607 mL	2.7213 mL
50 mM	0.0544 mL	0.2721 mL	0.5443 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481