


Cephapirin Sodium

Chemical Properties

CAS No. :	24356-60-3
Formula:	C ₁₇ H ₁₆ N ₃ O ₆ S ₂ .Na
Molecular Weight:	445.45
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Cefapirin sodium is the sodium salt of Cephapirin, an antimicrobial agent.
Targets(IC50)	Antibacterial, Antibiotic

Solubility Information

Solubility	DMSO: 20 mg/mL (44.9 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+90% Saline: 2 mg/mL (4.49 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2449 mL	11.2246 mL	22.4492 mL
5 mM	0.449 mL	2.2449 mL	4.4898 mL
10 mM	0.2245 mL	1.1225 mL	2.2449 mL
50 mM	0.0449 mL	0.2245 mL	0.449 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Tan JS, et al. Bacampicillin, Ampicillin, Cephalothin, and Cephapirin levels in human blood and interstitial fluid. Antimicrob Agents Chemother. 1979 Apr;15(4):510-2.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481