


5,7-dihydroxy-2-(4-[[trihydroxy-(hydroxymethyl)oxan-yl]oxy]phenyl)-dihydro-2H-1-benzopyran-4-one

Chemical Properties

CAS No. :

Formula: C₂₁H₂₂O₁₀

Molecular Weight: 434.4

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	5,7-dihydroxy-2-(4-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]phenyl)-3,4-dihydro-2H-1-benzopyran-4-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.302 mL	11.5101 mL	23.0203 mL
5 mM	0.4604 mL	2.302 mL	4.6041 mL
10 mM	0.2302 mL	1.151 mL	2.302 mL
50 mM	0.046 mL	0.2302 mL	0.4604 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481