

EN300-371838

## Chemical Properties

CAS No. : 2094748-09-9


Formula: C7H7ClN2O

Molecular Weight: 170.596

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

HCl



## Biological Description

## Description

EN300-371838, with CAS No. 2094748-09-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. EN300-371838 provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8617 mL	29.3083 mL	58.6166 mL
5 mM	1.1723 mL	5.8617 mL	11.7233 mL
10 mM	0.5862 mL	2.9308 mL	5.8617 mL
50 mM	0.1172 mL	0.5862 mL	1.1723 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481