


3-Methoxyfuran

Chemical Properties

CAS No. :	3420-57-3
Formula:	C ₅ H ₆ O ₂
Molecular Weight:	98.1
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Methoxyfuran shows Diels-Alder reactivity under forcing conditions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	10.1937 mL	50.9684 mL	101.9368 mL
5 mM	2.0387 mL	10.1937 mL	20.3874 mL
10 mM	1.0194 mL	5.0968 mL	10.1937 mL
50 mM	0.2039 mL	1.0194 mL	2.0387 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481