

Compound T130433

Chemical Properties


CAS No. :

Formula: C53H86O22

Molecular Weight: 1075.25

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2S,3R,4R,5R,6S)-2-[[[(2S,3R,4S,5S)-5-hydroxy-2-[[[(1S,2R,5R,7S,10R,14R,16S,20S)-16-hydroxy-2,6,6,10,16-pentamethyl-18-(2-methylprop-1-en-1-yl)-19,21-dioxahexacyclo[18.2.1.0 ^{1,14} .0 ^{2,11} .0 ^{5,10} .0 ^{15,20}]tricosan-7-yl]oxy]-4-[[[(2S,3R,4S,5S,6R)-3,4,5-trihydroxy-6-[[[(2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]methyl)oxan-2-yl]oxy]oxan-3-yl]oxy]-6-methyloxane-3,4,5-triol is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.930 mL	4.6501 mL	9.3002 mL
5 mM	0.186 mL	0.930 mL	1.860 mL
10 mM	0.093 mL	0.465 mL	0.930 mL
50 mM	0.0186 mL	0.093 mL	0.186 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481