

SARS-CoV-2-IN-20

Chemical Properties


CAS No. : 1260244-33-4

Formula: C₂₄H₃₀N₂O₂

Molecular Weight: 378.51

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	SARS-CoV-2-IN-20 (Compound 1a) is a highly effective SARS-CoV-2 inhibitor with an impressive EC ₅₀ value of 6.5 μM, demonstrating significant potential for advancing infectious disease research [1].
Targets(IC ₅₀)	Others,SARS-CoV

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6419 mL	13.2097 mL	26.4194 mL
5 mM	0.5284 mL	2.6419 mL	5.2839 mL
10 mM	0.2642 mL	1.321 mL	2.6419 mL
50 mM	0.0528 mL	0.2642 mL	0.5284 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481