


1H-Isoindol-1-one,2,3-dihydro-4-methoxy-

Chemical Properties

CAS No. : 366453-22-7
 Formula: C₉H₉NO₂
 Molecular Weight: 163.176
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	1H-Isoindol-1-one,2,3-dihydro-4-methoxy-, with CAS No. 366453-22-7, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 1H-Isoindol-1-one,2,3-dihydro-4-methoxy- provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.1282 mL	30.641 mL	61.282 mL
5 mM	1.2256 mL	6.1282 mL	12.2564 mL
10 mM	0.6128 mL	3.0641 mL	6.1282 mL
50 mM	0.1226 mL	0.6128 mL	1.2256 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481