


Vincanol

Chemical Properties

CAS No. :	19877-89-5
Formula:	C ₁₉ H ₂₄ N ₂ O
Molecular Weight:	296.41
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Vincanol is a potent voltage-gated Na ⁺ channels blocker.
Targets(IC50)	Others,Sodium Channel

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3737 mL	16.8685 mL	33.7371 mL
5 mM	0.6747 mL	3.3737 mL	6.7474 mL
10 mM	0.3374 mL	1.6869 mL	3.3737 mL
50 mM	0.0675 mL	0.3374 mL	0.6747 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Erdo SA, Molnár P, Lakics V, Bence JZ, Tömösközi Z. Vincamine and vincanol are potent blockers of voltage-gated Na⁺ channels. Eur J Pharmacol. 1996 Oct 24;314(1-2):69-73. PubMed PMID: 8957220.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481