


Clopidogrel carboxylic acid

Chemical Properties

CAS No. :	144457-28-3
Formula:	C ₁₅ H ₁₄ ClNO ₂ S
Molecular Weight:	307.79
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Clopidogrel carboxylic acid is a useful organic compound for research related to life sciences.
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.249 mL	16.2448 mL	32.4897 mL
5 mM	0.6498 mL	3.249 mL	6.4979 mL
10 mM	0.3249 mL	1.6245 mL	3.249 mL
50 mM	0.065 mL	0.3249 mL	0.6498 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481