

Mal-PEG4-VA

Chemical Properties

CAS No. : 1800456-31-8


Formula: C₂₆H₄₂N₄O₁₁

Molecular Weight: 586.63

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Mal-PEG4-VA, a noncleavable ADC linker featuring a Maleimide group, is utilized in the synthesis of antibody-drug conjugates.
Targets(IC50)	ADC Linker,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7047 mL	8.5233 mL	17.0465 mL
5 mM	0.3409 mL	1.7047 mL	3.4093 mL
10 mM	0.1705 mL	0.8523 mL	1.7047 mL
50 mM	0.0341 mL	0.1705 mL	0.3409 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481