

8,15-Isopimaradiene-1,7,14,18-tetrol, 18-Carboxyli

Chemical Properties


CAS No. :

Formula: C₂₇H₃₈O₈

Molecular Weight: 490.593

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	8,15-Isopimaradiene-1,7,14,18-tetrol, 18-Carboxyli is a useful organic compound for research related to life sciences and the catalog number is T130947.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0384 mL	10.1918 mL	20.3836 mL
5 mM	0.4077 mL	2.0384 mL	4.0767 mL
10 mM	0.2038 mL	1.0192 mL	2.0384 mL
50 mM	0.0408 mL	0.2038 mL	0.4077 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481