


Ro 15-2041

## Chemical Properties

CAS No. : 77378-53-1  
 Formula: C<sub>12</sub>H<sub>12</sub>BrN<sub>3</sub>O  
 Molecular Weight: 294.15  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ro 15-2041 is a selective platelet cAMP phosphodiesterase inhibitor with antithrombotic properties.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3996 mL	16.9981 mL	33.9963 mL
5 mM	0.6799 mL	3.3996 mL	6.7993 mL
10 mM	0.340 mL	1.6998 mL	3.3996 mL
50 mM	0.068 mL	0.340 mL	0.6799 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481