

## Compound T128334

## Chemical Properties

CAS No. :

Formula: C<sub>32</sub>H<sub>58</sub>O<sub>18</sub>

Molecular Weight: 730.8

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	5-[(5,7-dihydroxy-3-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]decanoyl)oxy]-3-[[3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy]decanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3684 mL	6.8418 mL	13.6836 mL
5 mM	0.2737 mL	1.3684 mL	2.7367 mL
10 mM	0.1368 mL	0.6842 mL	1.3684 mL
50 mM	0.0274 mL	0.1368 mL	0.2737 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481