

## trans-3,4-Methylenedioxcinnamyl alcohol

## Chemical Properties


CAS No. :

Formula: C<sub>10</sub>H<sub>10</sub>O<sub>3</sub>

Molecular Weight: 178.187

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	trans-3,4-Methylenedioxcinnamyl alcohol is a useful organic compound for research related to life sciences and the catalog number is T124617.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.612 mL	28.0599 mL	56.1199 mL
5 mM	1.1224 mL	5.612 mL	11.224 mL
10 mM	0.5612 mL	2.806 mL	5.612 mL
50 mM	0.1122 mL	0.5612 mL	1.1224 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481