


## Stemonidine

## Chemical Properties

CAS No. : 85700-47-6  
 Formula: C<sub>19</sub>H<sub>29</sub>NO<sub>5</sub>  
 Molecular Weight: 351.443  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Stemonidine is a natural product.
-------------	-----------------------------------

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8454 mL	14.2272 mL	28.4544 mL
5 mM	0.5691 mL	2.8454 mL	5.6909 mL
10 mM	0.2845 mL	1.4227 mL	2.8454 mL
50 mM	0.0569 mL	0.2845 mL	0.5691 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481