

## 2,3,4,5,6-pentahydroxyhexyl (2E,4E)-6-methylocta-2,4-dienoate

## Chemical Properties

CAS No. :

Formula: C<sub>15</sub>H<sub>26</sub>O<sub>7</sub>

Molecular Weight: 318.37

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	2,3,4,5,6-pentahydroxyhexyl (2E,4E)-6-methylocta-2,4-dienoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.141 mL	15.705 mL	31.410 mL
5 mM	0.6282 mL	3.141 mL	6.282 mL
10 mM	0.3141 mL	1.5705 mL	3.141 mL
50 mM	0.0628 mL	0.3141 mL	0.6282 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481