

Trepamotide

Chemical Properties

CAS No. : 1018833-44-7


Formula: C₅₈H₈₀N₁₀O₁₈

Molecular Weight: 1205.31

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

Trepamotide is a bioactive chemical.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8297 mL	4.1483 mL	8.2966 mL
5 mM	0.1659 mL	0.8297 mL	1.6593 mL
10 mM	0.083 mL	0.4148 mL	0.8297 mL
50 mM	0.0166 mL	0.083 mL	0.1659 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481