

HO-PEG-mal (MW 3400)

**Chemical Properties**

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	HO-PEG-mal (MW 3400) is a PEG-based linker for PROTACs, joining two essential ligands necessary for forming PROTAC molecules and enabling selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, composed of two distinct ligands connected by a linker—one for an E3 ubiquitin ligase and the other for the target protein—utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481