

## 1,2-Didehydrocryptotanshinone

## Chemical Properties

|                   |  |
|-------------------|--|
| CAS No. :         | 891854-92-5  |
| Formula:          | C19H18O3   |
| Molecular Weight: | 294.35   |
| Storage:          | Powder: -20°C for 3 years   In solvent: -80°C for 1 year<br><small>Actual storage temperature shall be subject to the COA.</small> |

## Biological Description

|               |   |
|---------------|---|
| Description   | 1,2-Didehydrocryptotanshinone is a useful organic compound for research related to life sciences and the catalog number is T123986. |
| Targets(IC50) | Others  |

## Preparing Stock Solutions

|       | 1mg       | 5mg        | 10mg       |
|-------|-----------|------------|------------|
| 1 mM  | 3.3973 mL | 16.9866 mL | 33.9732 mL |
| 5 mM  | 0.6795 mL | 3.3973 mL  | 6.7946 mL  |
| 10 mM | 0.3397 mL | 1.6987 mL  | 3.3973 mL  |
| 50 mM | 0.0679 mL | 0.3397 mL  | 0.6795 mL  |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481