


3'-Hydroxydehydroaglaiastatin

Chemical Properties

CAS No. :	259143-58-3
Formula:	C ₃₁ H ₂₈ N ₂ O ₇
Molecular Weight:	540.56
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3'-Hydroxydehydroaglaiastatin is a natural product for research related to life sciences. The catalog number is TN2948 and the CAS number is 259143-58-3.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8499 mL	9.2497 mL	18.4993 mL
5 mM	0.370 mL	1.8499 mL	3.6999 mL
10 mM	0.185 mL	0.925 mL	1.8499 mL
50 mM	0.037 mL	0.185 mL	0.370 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481