

Compound N038-0225

Chemical Properties


CAS No. :

Formula: C₂₃H₃₂O₄

Molecular Weight: 372.505

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N038-0225 is a useful organic compound for research related to life sciences and the catalog number is T131279.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6845 mL	13.4225 mL	26.8449 mL
5 mM	0.5369 mL	2.6845 mL	5.369 mL
10 mM	0.2684 mL	1.3422 mL	2.6845 mL
50 mM	0.0537 mL	0.2684 mL	0.5369 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481