


[4-(Tetrahydropyran-4-yloxy)phenyl]methanol

Chemical Properties

CAS No. :	892501-95-0
Formula:	C ₁₂ H ₁₆ O ₃
Molecular Weight:	208.26
Storage:	Store under nitrogen Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	[4-(Tetrahydropyran-4-yloxy)phenyl]methanol, with CAS No. 892501-95-0, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. [4-(Tetrahydropyran-4-yloxy)phenyl]methanol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8017 mL	24.0085 mL	48.0169 mL
5 mM	0.9603 mL	4.8017 mL	9.6034 mL
10 mM	0.4802 mL	2.4008 mL	4.8017 mL
50 mM	0.096 mL	0.4802 mL	0.9603 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481