

cIAP1 Ligand-Linker Conjugates 9

Chemical Properties

CAS No. :


Formula: C37H48N4O7

Molecular Weight: 660.8

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	cIAP1 Ligand-Linker Conjugates 9 feature an IAP ligand specifically targeting the E3 ubiquitin ligase, along with a PROTAC linker. It is valuable tools for developing SNIPERs [1].
Targets(IC50)	Others,E3 Ligase Ligand-Linker Conjugates,PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5133 mL	7.5666 mL	15.1332 mL
5 mM	0.3027 mL	1.5133 mL	3.0266 mL
10 mM	0.1513 mL	0.7567 mL	1.5133 mL
50 mM	0.0303 mL	0.1513 mL	0.3027 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Itoh Y, et al. Design, synthesis and biological evaluation of nuclear receptor-degradation inducers. Bioorg Med Chem. 2011 Nov 15;19(22):6768-78.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481