


5-Methyl-2-thiouridine

Chemical Properties

CAS No. :	32738-09-3
Formula:	C10H14N2O5S
Molecular Weight:	274.29
Storage:	Store at low temperature Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	5-Methyl-2-thiouridine is a purine nucleoside analog, applicable for biochemical experiments and drug synthesis research.
Targets(IC50)	Nucleoside Antimetabolite/Analog

Solubility Information

Solubility	DMSO: 200 mg/mL (729.16 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6458 mL	18.2289 mL	36.4578 mL
5 mM	0.7292 mL	3.6458 mL	7.2916 mL
10 mM	0.3646 mL	1.8229 mL	3.6458 mL
50 mM	0.0729 mL	0.3646 mL	0.7292 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481