


## Soyasaponin A1

## Chemical Properties

CAS No. :	78693-94-4
Formula:	C <sub>59</sub> H <sub>96</sub> O <sub>29</sub>
Molecular Weight:	1269.39
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Soyasaponin A1 is a useful organic compound for research related to life sciences. The catalog number is T124217 and the CAS number is 78693-94-4.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7878 mL	3.9389 mL	7.8778 mL
5 mM	0.1576 mL	0.7878 mL	1.5756 mL
10 mM	0.0788 mL	0.3939 mL	0.7878 mL
50 mM	0.0158 mL	0.0788 mL	0.1576 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481