


Symphonone I

Chemical Properties

CAS No. :	1235774-18-1
Formula:	C ₃₈ H ₄₈ O ₆
Molecular Weight:	600.796
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Symphonone I is a natural product for research related to life sciences. The catalog number is TN5908 and the CAS number is 1235774-18-1.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6644 mL	8.3222 mL	16.6445 mL
5 mM	0.3329 mL	1.6644 mL	3.3289 mL
10 mM	0.1664 mL	0.8322 mL	1.6644 mL
50 mM	0.0333 mL	0.1664 mL	0.3329 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481