


BN80933

Chemical Properties

CAS No. : 214348-10-4
 Formula: C₂₉H₃₄N₄O₃S
 Molecular Weight: 518.67
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	BN80933 is a dual inhibitor of lipid peroxidation and neuronal nitric oxide synthase; Neuroprotectant.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.928 mL	9.640 mL	19.2801 mL
5 mM	0.3856 mL	1.928 mL	3.856 mL
10 mM	0.1928 mL	0.964 mL	1.928 mL
50 mM	0.0386 mL	0.1928 mL	0.3856 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481