


(S)-GNA-A(Bz)-phosphoramidite

Chemical Properties

CAS No. : 851050-24-3
 Formula: C₄₅H₅₀N₇O₆P
 Molecular Weight: 815.9
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-GNA-A(Bz)-phosphoramidite is a Nucleoside Phosphoramidite.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2256 mL	6.1282 mL	12.2564 mL
5 mM	0.2451 mL	1.2256 mL	2.4513 mL
10 mM	0.1226 mL	0.6128 mL	1.2256 mL
50 mM	0.0245 mL	0.1226 mL	0.2451 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481