


## 2,16-Kauranediol

## Chemical Properties

CAS No. :	34302-37-9
Formula:	C <sub>20</sub> H <sub>34</sub> O <sub>2</sub>
Molecular Weight:	306.49
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2,16-Kauranediol is a natural product for research related to life sciences. The catalog number is TN2676 and the CAS number is 34302-37-9.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2627 mL	16.3137 mL	32.6275 mL
5 mM	0.6525 mL	3.2627 mL	6.5255 mL
10 mM	0.3263 mL	1.6314 mL	3.2627 mL
50 mM	0.0653 mL	0.3263 mL	0.6525 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481