


WEHL-03

## Chemical Properties

CAS No. : 15827-80-2  
 Formula: C<sub>6</sub>H<sub>9</sub>N<sub>3</sub>  
 Molecular Weight: 123.16  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	WEHL-03 (4-AMINO-3-CYANO-1,2,5,6-TETRAHYDROPYRIDINE) is an active chemical.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: 60 mg/mL (487.17 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.1195 mL	40.5976 mL	81.1952 mL
5 mM	1.6239 mL	8.1195 mL	16.239 mL
10 mM	0.812 mL	4.0598 mL	8.1195 mL
50 mM	0.1624 mL	0.812 mL	1.6239 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481