


## [Des-Arg9]-Bradykinin

## Chemical Properties

CAS No. :	15958-92-6
Formula:	C44H61N11O10
Molecular Weight:	904.02
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	[Des-Arg9]-Bradykinin ([Des-Arg9]-Bradykinin(2TFA))(2TFA) is an agonist of Bradykinin (B1) receptor.
Targets(IC50)	Bradykinin Receptor

## Solubility Information

Solubility	DMSO: 10 mM, Sonication is recommended. ( $< 1$ mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1062 mL	5.5309 mL	11.0617 mL
5 mM	0.2212 mL	1.1062 mL	2.2123 mL
10 mM	0.1106 mL	0.5531 mL	1.1062 mL
50 mM	0.0221 mL	0.1106 mL	0.2212 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Leeb-Lundberg, L. M F . International Union of Pharmacology. XLV. Classification of the Kinin Receptor Family: from Molecular Mechanisms to Pathophysiological Consequences[J]. Pharmacological Reviews, 2005, 57(1):27-77.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481