

## HDAC1-IN-3

## Chemical Properties

CAS No. :

Formula: C<sub>22</sub>H<sub>24</sub>ClN<sub>7</sub>O<sub>2</sub>

Molecular Weight: 453.92

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	HDAC1-IN-3 is a potent inhibitor of Pf HDAC1, exhibiting a significant in vivo killing effect on all life cycles of the parasite and displaying antimalarial activity in both wild-type and multidrug-resistant strains of the parasite.
Targets(IC50)	Others,HDAC,Parasite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.203 mL	11.0152 mL	22.0303 mL
5 mM	0.4406 mL	2.203 mL	4.4061 mL
10 mM	0.2203 mL	1.1015 mL	2.203 mL
50 mM	0.0441 mL	0.2203 mL	0.4406 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481