


NU-8165

Chemical Properties

CAS No. : 849723-20-2
 Formula: C₂₄H₂₂ClNO₃
 Molecular Weight: 407.89
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	NU-8165 is an MDM2-p53 protein-protein interaction inhibitor.
Targets(IC50)	Others,MDM-2/p53

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4516 mL	12.2582 mL	24.5164 mL
5 mM	0.4903 mL	2.4516 mL	4.9033 mL
10 mM	0.2452 mL	1.2258 mL	2.4516 mL
50 mM	0.049 mL	0.2452 mL	0.4903 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481