


## 2-Amino-1-beta-D-arabinofuranosyl-5-methyl-4(1H)-pyrimidinone

## Chemical Properties

CAS No. : 10212-31-4  
 Formula: C<sub>10</sub>H<sub>15</sub>N<sub>3</sub>O<sub>5</sub>  
 Molecular Weight: 257.24  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	Nucleoside Derivatives - Arabino-nucleosides, Iso-cytidine derivatives
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8874 mL	19.4371 mL	38.8742 mL
5 mM	0.7775 mL	3.8874 mL	7.7748 mL
10 mM	0.3887 mL	1.9437 mL	3.8874 mL
50 mM	0.0777 mL	0.3887 mL	0.7775 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481