

JNJ-40264796

Chemical Properties


CAS No. :

Formula: C₂₂H₂₁N₃O₂

Molecular Weight: 359.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	JNJ-40264796 is a negative control for JNJ-40068782.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7822 mL	13.9109 mL	27.8218 mL
5 mM	0.5564 mL	2.7822 mL	5.5644 mL
10 mM	0.2782 mL	1.3911 mL	2.7822 mL
50 mM	0.0556 mL	0.2782 mL	0.5564 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481