

## Digitoxin

## Chemical Properties

CAS No. : 71-63-6

Formula: C<sub>41</sub>H<sub>64</sub>O<sub>13</sub>


Molecular Weight: 764.94

Storage:

Keep away from direct sunlight, Keep away from moisture, Store at low temperature

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Digitoxin (Digitoxinum) is an effective Na <sup>+</sup> /K <sup>+</sup> -ATPase inhibitor (EC <sub>50</sub> value: 0.78 μM).
Targets(IC <sub>50</sub> )	Apoptosis,ATPase,Bcl-2 Family,Calcium Channel,Caspase,HSV

## Solubility Information

Solubility	DMSO: 100 mg/mL (130.73 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 4 mg/mL (5.23 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3073 mL	6.5365 mL	13.0729 mL
5 mM	0.2615 mL	1.3073 mL	2.6146 mL
10 mM	0.1307 mL	0.6536 mL	1.3073 mL
50 mM	0.0261 mL	0.1307 mL	0.2615 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Bossuyt X, et al. J Pharmacol Exp Ther. 1996 Mar;276(3):891-6.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481