

Acetyl Tetrapeptide-5 acetate

Chemical Properties

CAS No. :


Formula: C₂₂H₃₂N₈O₉

Molecular Weight: 552.55

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Acetyl Tetrapeptide-5 acetate helps reduce eye puffiness, improve skin elasticity as well as overall smoothness.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 2.5 mM, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8098 mL	9.049 mL	18.0979 mL
5 mM	0.362 mL	1.8098 mL	3.6196 mL
10 mM	0.181 mL	0.9049 mL	1.8098 mL
50 mM	0.0362 mL	0.181 mL	0.362 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481