

Ibuprofen Lysine

Chemical Properties


CAS No. : 57469-77-9

Formula: C₁₉H₃₂N₂O₄

Molecular Weight: 352.47

Storage: Store at low temperature, Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ibuprofen Lysine (Neoprofen) is a non-steroidal anti-inflammatory drug.
Targets(IC50)	Apoptosis, Parasite, COX

Solubility Information

Solubility	H ₂ O: 250 mg/mL (709.28 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8371 mL	14.1856 mL	28.3712 mL
5 mM	0.5674 mL	2.8371 mL	5.6742 mL
10 mM	0.2837 mL	1.4186 mL	2.8371 mL
50 mM	0.0567 mL	0.2837 mL	0.5674 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Ferrero-Cafiero JM, et al. Int J Clin Pharmacol Ther. 2015 Nov;53(11):972-9.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481