


Thevetin A

Chemical Properties

CAS No. :	37933-66-7
Formula:	C42H64O19
Molecular Weight:	872.955
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Thevetin A is a useful organic compound for research related to life sciences. The catalog number is T124930 and the CAS number is 37933-66-7.
Targets(IC50)	Others,Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1455 mL	5.7276 mL	11.4553 mL
5 mM	0.2291 mL	1.1455 mL	2.2911 mL
10 mM	0.1146 mL	0.5728 mL	1.1455 mL
50 mM	0.0229 mL	0.1146 mL	0.2291 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481