


Phthalimide, N-(5-(p-aminophenoxy)pentyl)-3-nitro-

Chemical Properties

CAS No. :	108368-14-5
Formula:	C19H19N3O5
Molecular Weight:	369.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phthalimide, N-(5-(p-aminophenoxy)pentyl)-3-nitro- is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7073 mL	13.5366 mL	27.0731 mL
5 mM	0.5415 mL	2.7073 mL	5.4146 mL
10 mM	0.2707 mL	1.3537 mL	2.7073 mL
50 mM	0.0541 mL	0.2707 mL	0.5415 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481