


## IS-741 potassium

## Chemical Properties

CAS No. :	141284-74-4
Formula:	C <sub>15</sub> H <sub>19</sub> F <sub>3</sub> KN <sub>3</sub> O <sub>3</sub> S
Molecular Weight:	417.49
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	IS-741 potassium is a phospholipase A2 inhibitor.
Targets(IC50)	Others,Integrin,Phospholipase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3953 mL	11.9763 mL	23.9527 mL
5 mM	0.4791 mL	2.3953 mL	4.7905 mL
10 mM	0.2395 mL	1.1976 mL	2.3953 mL
50 mM	0.0479 mL	0.2395 mL	0.4791 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481