


## Biotin-PEG1-azide

## Chemical Properties

CAS No. :	1204085-48-2
Formula:	C <sub>14</sub> H <sub>24</sub> N <sub>6</sub> O <sub>3</sub> S
Molecular Weight:	356.44
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Biotin-PEG1-azide is a 1-unit polyethylene glycol (PEG) linker used in the synthesis of antibody-drug conjugates (ADCs).
Targets(IC50)	ADC Linker
In vitro	ADCs consist of an antibody linked to an ADC cytotoxin via an ADC linker[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8055 mL	14.0276 mL	28.0552 mL
5 mM	0.5611 mL	2.8055 mL	5.611 mL
10 mM	0.2806 mL	1.4028 mL	2.8055 mL
50 mM	0.0561 mL	0.2806 mL	0.5611 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Beck A, et al. Strategies and challenges for the next generation of antibody-drug conjugates. Nat Rev Drug Discov. 2017 May;16(5):315-337.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481