

CG 428-Neg

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Negative control for CG 428.
Targets(IC50)	Others

Reference

Chen et al (2020) Discovery of first-in-class potent and selective tropomyosin receptor kinase degraders. J.Med. Chem. 63 14562 PMID: 33058680

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481